



LEARN MORE!



BELPYCNO

AUTOMATIC HIGH-PRECISION PYCNOMETER

**DENSITY MEASUREMENT OF 3D PRINTER
ADDITIVE MOLDING MATERIALS
BY GAS DISPLACEMENT METHOD**

Highly accurate, quick, and easy to use

Density measurement of 3D printer additive molding material by gas displacement method

Evaluation of true density by gas displacement method, which is quick and easy to measure, is effective for understanding the strength of layered objects with complex shapes created by 3D printers



3D PRINTER LAYERED OBJECTS



LOW DENSITY ... HIGH POROSITY ... FRAGILE STRUCTURE
HIGH DENSITY ... LOW POROSITY ... ROBUST CONSTRUCTION



Understanding density is essential to structure determination



EXAMPLE: SIMPLE LAYERED MOLDING BODY



EXAMPLE: COMPLEX LAYERED MOLDINGS



Stacked molding
body image
by 3D printing



FIRING (DEBINDING) AND SINTERING

ELECTRIC FURNACE: CARBORITE



SIMPLE SHAPE AND FEW VOIDS

| Density (P) = W (Weight) / V (Volume)

| V (Volume) = L (Length) x W (Width) x H (Height)




REGARDLESS OF COMPLEX SHAPES AND LARGE OR SMALL VOIDS

| Density evaluation by Archimedes method (displacement in water)

Problem

Liquid must be sufficiently wetted to the voids of the molding body (long stabilization time)



| Density evaluation by gas displacement method

Advantage

Short measurement time with or without voids

BELPYCNO SERIE



BELPYCNO
Dedicated to metal powders
and other powders



BELPYCNO L
For density evaluation of
large-volume moldings



SPECIFICATIONS

Measurement principle	Gas displacement method
Sample cell volume	10 cm ³ , 3.5 cm ³ , 1 cm ³
Measurement accuracy	0.03% F.S + 0.03% R
Repeatability	0.02% F.S
Pretreatment	Gas purge, Flow, Vacuum (option)
Pretreatment pressure	0~145 kPa (Gauge)
Measurement pressure	145 kPa (Gauge)
Measurement repetition number of time	Max. 100 times
Mean number of times	Max. 100 times
Measurement temperature	Room temperature (Circulator (Option)) : 15 to 50 °C
Calibration method	Automatic calibration with calibration sphere
Interface Input	RS232C (communication with a electric balance)
Utility gas	He, inert gases: pres. 1.5bar (Gauge), 1/8" Swagelok joint
Allowed gases	He, N ₂ , inert gases
Power supply	AC 90~250 V / 200 W
Dimensions, Weight	270 (W) × 170 (H) × 300 (D) mm, 8 kg
Recommended Monitor	Full HD monitors

DATA SHEET

Measurement result output is in text-file format and can be printed in report form. Measurement data can also be edited by Microsoft Excel.

▼Detail of measurement result



▼List of measurement results



- Specifications and appearance of the products listed are subject to change without notice.
- Products (goods and services) described in the catalog, depending on the destination and application, might be applicable to export regulations, etc. by the "Foreign Exchange and Foreign Trade Control Law".
- In response to the review of the Japanese government regarding the export of products (goods and services), permission and approval, and the like, must be obtained according to the regulations.

MICROTRAC

PARTICLE CHARACTERIZATION


Microtrac Inc.
11 Penns Trail • Newtown, PA 18940 • USA
Phone: +1 866 473 8724 • marketing@microtrac.com

Microtrac Retsch GmbH
Retsch-Allee 1-5 • 42781 Haan • Germany
Phone: +49 2104 2333 300 • info@microtrac.com

MicrotracBEL Corp.
8-2-52 Nanko Higashi, Suminoe-ku • Osaka 559-0031 • Japan
Phone: +81 6-6655-0362 • international@microtrac-bel.com

Microtrac Formulaction SAS
3-5 rue Paule Raymondis • 31200 Toulouse • France
Phone: +33 (0)5 62 89 29 29 • contact.fr@mtf.verder.com

www.microtrac.



part of **VERDER**